

## PROJECT TWO: MILESTONE 1 – COVER PAGE

Team Number: 

THURS- 12
--------------

Please list full names and MacID's of all *present* Team Members

Full Name:	MacID:
Yuvraj Sandhu	Sandhuy
Eric Hitsman	hitsmane
Chengyao Liu	liuc169
Couper Smith	smitc25

## MILESTONE 1 (STAGE 1) – PRE-PROJECT ASSIGNMENT

Team Number: 

Thurs-12
----------

You should have already completed this task individually prior to Design Studio 7.

1. Copy-and-paste each team member's list of objectives, constraints and functions on the following pages (1 team member per page)
  - a. Be sure to indicate each team member's Name and MacID

We are asking that you submit your work on both worksheets. It does seem redundant, but there are valid reasons for this:

- Each team member needs to submit their list of objectives, constraints and functions with the **Milestone One Individual Worksheets** document so that it can be **graded**
- Compiling your individual work into this **Milestone One Team Worksheets** document allows you to readily access your team member's work
  - This will be especially helpful when completing **Stage 2** of the milestone

Team Number: Thurs-12

Name:Eric Hitsman	MacID:hitsmane
<i>Objectives</i> <ul style="list-style-type: none"><li>• Easy for machine to access</li><li>• Holds tools securely</li><li>• Tools return to same spot every time</li></ul> <i>Constraints</i> <ul style="list-style-type: none"><li>• Must be secure</li><li>• Must not exceed 350g</li><li>• Container must be smaller than grabber</li></ul> <i>Functions</i> <ul style="list-style-type: none"><li>• Hold surgical tools</li><li>• Sort containers</li></ul>	

Team Number: Thurs-12

Name: Chengyao	MacID: liuc169
<i>Objectives</i> <ul style="list-style-type: none"><li>Easy to be picked up by the robotic arm.</li><li>Able to hold the surgical tools securely.</li></ul> <i>Constraints</i> <ul style="list-style-type: none"><li><math>M \leq 350g</math></li></ul> <i>Functions</i> <ul style="list-style-type: none"><li>Able to securely hold surgical instrument.</li><li>Held by robotic arm stably.</li><li>Classify the surgical tools by colors.</li></ul>	

Team Number: 

Thurs-12
----------

Name: Yuvraj Sandhu	MacID: Sandhuy
<i>Objectives</i> <ul style="list-style-type: none"><li>• Easily be able to carry</li><li>• Fits within the claw's dimensions</li></ul> <i>Constraints</i> <ul style="list-style-type: none"><li>• Safe to utilize</li><li>• Must be able to carry up to 1 kilogram</li><li>• Cannot weigh larger than 350g</li><li>• Holds surgical tool in place without moving</li></ul> <i>Functions</i> <ul style="list-style-type: none"><li>• Firmly holds objects mid-air</li><li>• Can store solids and liquids without any leakage</li><li>• Transports objects to correct</li><li>• Allows steam to sanitize surgical tools</li></ul>	

Team Number: Thurs-12

Name: Couper Smith	MacID: smitc25
<b>Containers:</b> <i>Objectives</i> <ul style="list-style-type: none"><li>• Made of material with high coefficient of friction (to prevent slipping)</li><li>• Sturdy design (can survive drops and not spill contents)</li></ul> <i>Constraints</i> <ul style="list-style-type: none"><li>• Must NOT have a length or width greater than 150mm (proximal end limits of end-effector)</li><li>• Size of the container must fit inside the autoclave (different dimensions depending on assigned scenario)</li><li>• All features must be greater than 4mm (minimum dimension to be properly fabricated)</li><li>• Must not exceed 350g</li></ul> <i>Functions</i> <ul style="list-style-type: none"><li>• Able to securely hold the surgical tool in place during transfer</li><li>• Must allow sterilization to occur (ie steam must be able to penetrate the container)</li></ul>	<b>Robotic Arm/Autoclave Bins:</b> <i>Objectives</i> <ul style="list-style-type: none"><li>• Reliable (ie is able to ID the container and put it in the right location each time)</li><li>• Require minimal human intervention Must allow sterilization to occur (ie steam must be able to penetrate the container)</li></ul> <i>Constraints</i> <ul style="list-style-type: none"><li>• Must provide enough force to move the container but not enough force to break the container</li></ul> <i>Functions</i> <ul style="list-style-type: none"><li>• Able to identify the type of container placed on the pick-up platform</li><li>• Able to pick up the containers and securely move them to the correct container</li></ul>

## MILESTONE 1 (STAGE 2) – LIST OF OBJECTIVES, CONSTRAINTS, AND FUNCTIONS

Team Number: Thurs-12

1. As a team, create a final a list of objectives, constraints, and functions in the table below.

- Use your individual *Pre-Project Assignment* to build your team's final list
- The exact number you should have depends on what information you have gathered from the Project Pack.

Objectives	Constraints	Functions
Made of material with high coefficient of friction (to prevent slipping)	Must not exceed 350g	Must allow sterilization to occur (ie steam must be able to penetrate the container)
Sturdy design (can survive drops and not spill contents)	Must NOT have a length or width greater than 150mm (proximal end limits of end-effector)	Securely hold the surgical tool in place during transfer
Easily able to be carried	All features must be greater than 4mm (minimum dimension to be properly fabricated)	Classify the surgical tools by colors and size
Holds tools securely	Size of the container must fit inside the autoclave (different dimensions depending on assigned scenario)	

2. What is the primary function of the entire system?

Allow Transport of Surgical Tools

3. What are the secondary functions?

Sterilizes Tools
Correctly Identifies Container and End Location
Securely hold the surgical devices in place



## MILESTONE 1 (STAGE 3) – MORPHOLOGICAL ANALYSIS

Team Number: **Thurs-12**

1. Identify multiple means to perform the secondary functions that your team came up with during Stage 1 of this milestone. One sub-function (pick up) is already listed for you. The other two sub-functions are for your team to choose.

→ Make sure that every mean for the “pick-up” sub-function assumes that the end effector of the robot arm is a gripper. The means for your other sub-functions do not need to follow this assumption.

Function	Means					
Pick up	Grasp (from the sides)	Lift from bottom	Hoop on container	Suction Cups	Magnetic	Handle
Securely holds devices	Container	Indents in Container	High Friction Claw Material	Magnetic	Clips	Velcro
Correctly identifies end location	Camera	Coordinates system	Colour sensor	Depth sensor	Pre-Planned Code	

## MILESTONE 1 (STAGE 4) – CONCEPT SKETCHES

Team Number: 

Thurs-12
----------

Complete this worksheet *after* having completed stage 3 as a team **and** after having **individually** created your concept sketches.

1. Each team member should copy-and-paste the photo of their individual concept sketches in the space indicated on the following pages
  - The photo's should be the same one you included in the **Milestone One Individual Worksheets** document
  - Be sure to include your **Team Number** on each page
  - Be sure each team member's **Name** and **MacID** are included with each sketch

We are asking that you submit your work on both worksheets. It does seem redundant, but there are valid reasons for this:

- Each team member needs to submit their sketch with the **Milestone One Individual Worksheets** document so that it can be **graded**
- Compiling your individual work into this **Milestone One Team Worksheets** document allows you to readily access your team member's work

Name: Yuvraj Sandhu


MacID: Sandhuy

Yuvraj Sandhu  
Sandhuy  
THUR-12


Design Features:

- 2 Clips for security
- Suctions cups to firmly grasp the indented cups on the container

TOP VIEW




Multi-view




2 cup indents to maximize grip of the suction cup since the indent increases the suction  
2 clips on both sides to lock the claw into place.

LEFT VIEW



claw




suction cups

Yuvraj Sandhu  
Sandhuy  
THURS-12


Design Features:

- 2 Gripping handles for gripping in multiple angles.
- Indents for better grippage of the claw.
- 3 secure positions for transportation
- Rubber claw pads to increase grip on container.


TOP VIEW




Multi-VIEW




FRONT VIEW



SIDE VIEW



Claw Design





Rubber pads on claw to maximize grip

Name: Chengyao Liu

MacID: liuc169

①




front


dig 2 holes to make clipping stably.

Chengyao Liu  
liuc169  
Thurs-12.

②



suction Cup



left

put the  
tool at here,  
the sloop can make  
sure the tool won't slide out.

this layer  
is pretty smooth,  
polish or take a  
piece of glass on it.

Chengyao Liu  
liuc169  
Thurs-12


Team Number: Thurs-12

Name: Eric Hitsman

MacID: hitsmane

①

Indents



multiview




Front View


Eric Hitsman  
Thurs-12  
hitsmane

②

Magnets



multiview




Front view

Eric Hitsman  
Thurs-12  
hitsmane


Team Number: Thurs-12

Name: Couper Smith

MacID: smitc25



Thurs-12  
Couper Smith  
smitc25



Thurs-12  
Couper Smith  
smitc25